

Recombinant MTB PDZ Protein

Background

Recombinant MTB PDZ Protein is produced in *E.coli*.

Product Details

Product Code	RMTB004
Molecular weight	8 kDa
Sequence	AFL GLGVVDNNGN GARVQRVVGSAAPAASLGISG DVITAVDVGAP INSATAMADALNGHHPGDVIVSITWQTKSGG TRTGNTVLAEGPPA
Host	E.Coli
Purity/ Purification	> 96% Purity, purified by Ni NTA and determined by SDS-PAGE.
Form	Lyophilized powder. Dilute to 1 mg/mL in PBS
Purification Tag	C- terminus 6 X His tag
Applications	Tuberculosis research studies
Preparation and Storage	Store at -20 degree C, for extended storage
ISO Certification	Manufactured in an ISO 9001:20015 Certified Laboratory.

SDS-PAGE MTB PDZ Protein



Ordering Information:
GeNext Genomics Pvt. Ltd.
103, Abhyankar Nagar,
Nagpur
Ph.No.08888803973, 08380084087