

## MTB PepQ Monoclonal Antibody


### Background

Anti MTB PepQ MAb-MAB3H6D10 is monoclonal antibody which recognizes Mycobacterium tuberculosis PepQ protein. The MTB PepQ protein is essential for the survival and pathogenesis of the organism. This monoclonal antibody is generating using mouse as a host.

### Product Details

Catalog number	MAB3H6D10
Quantity	100µg
Size	150kda
Host	Mouse
Clone	6D10
Clonality	Monoclonal
Immunogen	Recombinant MTB PepQ Protein
Conjugate / Tag	Unconjugated
Purity	Protein A purified
Form	Lyophilized
Dilution Buffer (to be used)	1X PBS as per the concentration required for assay.
Antibody product type	Primary Antibody
Applications	ELISA (1µg/ml), WB (2µg/ml)

### Western Blot image of MTB PepQ Antibody



#### Storage:

For long term storage store at -20°C in small aliquots at a concentration of 0.1mg/ml to 0.5mg/ml to prevent freeze-thaw cycles.