

MTB PepA Monoclonal Antibody


Background

Anti MTB PepA MAb-MAB51DG9 is monoclonal antibody which recognizes Mycobacterium tuberculosis PepA protein. The MTB PepA protein is essential for the survival and pathogenesis of the organism. This monoclonal antibody is generating using mouse as a host.

Product Details

Catalog number	MAB51DG9
Quantity	100µg
Size	150kda
Host	Mouse
Clone	G9
Clonality	Monoclonal
Immunogen	Recombinant MTB PepA Protein
Conjugate / Tag	Unconjugated
Purity	Protein A purified
Form	Lyophilized
Dilution Buffer (to be used)	1X PBS as per the concentration required for assay.
Antibody product type	Primary Antibody
Applications	ELISA (1µg/ml), WB (2µg/ml)

Western Blot image of MTB PepA Antibody



Storage:

For long term storage store at -20°C in small aliquots at a concentration of 0.1mg/ml to 0.5mg/ml to prevent freeze-thaw cycles.