

Recombinant MTB PDZ Protein

Product Details

Catalog no.	RMTB004
Molecular weight	8kDa
Uniprot ID	UPI000E365BEC
Amino acid sequence	A1-A87
Host	E.coli
Formulation	50mM Tris pH 8.0 + 300mM NaCl + 1% Glycerol
Conjugate/ Tag	C-terminus 6 × His tag
Purity	> 90% Purity, purified by Ni NTA and determined by SDS-PAGE
Form	Liquid/ Lyophilized
Storage	Store at -20°C and avoid freeze thaw. Stable for 12 months from the date of receipt if kept at recommended temperature
Application	Tuberculosis research
Suggested working dilution	The optimal concentration should be determined for each specific application

SDS PAGE Gel image of MTB DapE protein

