

Monoclonal Anti STAT Antibody

The antibody is prepared against Recombinant fragment, corresponding to a region within amino acids 529-747 of Human STAT3. STAT3 is a member of the STAT Protein family. In response to cytokines and growth factors, STAT3 is phosphorylated by receptor-associated Janus Kinase (JAK), form homo- or heterodimers, and translocate to the cell nucleus where they act as transcription activator.

- **Product code:** MABSTATA10
- **Qty:**
- **Lot No.:**
- **Expiry:**

Package Contents:

Product Specifications

Description: Purified Mouse Immunoglobulin Clone 1A10

Specificity: It binds to N-Terminal domain of human STAT3.

Purification: Purified through Protein A Column

Clonality: Monoclonal

Target: Human STAT 3

Applications: ELISA, IHC, WB

