

Recombinant TNF α Antibody

Product Details

Catalog no.	RTNVH_002
Quantity	100 μ g
Size	15kDa
Clonality	Monoclonal
Subtype	VHH
Immunogen	Recombinant Human TNF α protein
Conjugate/ Tag	Unconjugated
Purification	Ni-NTA Purified
Buffer	Phosphate buffered saline (pH 7.4)
Form	Lyophilized/ Liquid
Resuspend	1X PBS or MilliQ Water
Application	ELISA, WB
Usage recommendations	ELISA- 1 μ g/ml, WB-1 μ g/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

WB showing reactivity of Anti-TNF- α VHH antibody with 17 kDa fraction of TNF- α protein

