



Sphingomyelinase haemolysins Monoclonal Antibody

Product Details

Catalog no.	MABCH310
Quantity	100µg
Size	150kDa
Host	Mouse
Clone	3H10
Clonality	Monoclonal
Subtype	IgG1
Immunogen	Recombinant SPH2 protein
Conjugate/ Tag	Unconjugated
Purification	Purified by Protein G affinity chromatography
Buffer	Phosphate buffered saline (pH 7.4)
Form	Lyophilized/ Liquid
Resuspend	Resuspend in MilliQ water
Application	ELISA, WB
Usage recommendations	ELISA- 1µg/ml, WB-1µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

WB image of Sphingomyelinase haemolysins Monoclonal Antibody (MABCH3C10)

