



CP4 EPSPS Polyclonal Antibody

Product Details

Catalog no.	PABHTBT02
Quantity	100µg
Size	150kDa
Host	Rabbit
Clonality	Polyclonal
Titre	30ng/ml
Immunogen	Recombinant CP4 EPSPS protein
Conjugate/ Tag	Unconjugated
Purification	Purified by ion exchange chromatography
Buffer	Phosphate buffered saline (pH 7.4)
Form	Lyophilized/ Liquid
Resuspend	Resuspend in MilliQ water
Application	ELISA
Usage recommendations	ELISA- 1µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

CP4-EPSPS Polyclonal Antibody Titre & Sandwich ELISA

