

IL 6 Receptor Monoclonal Antibody

Product Details

Catalog no.	MAB1F3
Quantity	100µg
Size	150kDa
Host	Mouse
Clone	1F3
Clonality	Monoclonal
Subtype	IgG2a
Immunogen	Recombinant Human IL6 receptor protein
Conjugate/ Tag	Unconjugated
Purification	Purified by protein A affinity chromatography
Buffer	Phosphate buffered saline (pH 7.4)
Form	Lyophilized/ Liquid
Resuspend	Resuspend in MilliQ water
Application	ELISA, WB
Usage recommendations	ELISA- 1µg/ml, WB-1µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

WB image of IL 6 Receptor Monoclonal Antibody (MAB1F3)

