

## IL 6 Receptor Monoclonal Antibody

### Product Details

<b>Catalog no.</b>	MAB1F3
<b>Quantity</b>	100µg
<b>Size</b>	150kDa
<b>Host</b>	Mouse
<b>Clone</b>	1F3
<b>Clonality</b>	Monoclonal
<b>Subtype</b>	IgG2a
<b>Immunogen</b>	Recombinant Human IL6 receptor protein
<b>Conjugate/ Tag</b>	Unconjugated
<b>Purification</b>	Purified by protein A affinity chromatography
<b>Buffer</b>	Phosphate buffered saline (pH 7.4)
<b>Form</b>	Lyophilized/ Liquid
<b>Resuspend</b>	Resuspend in MilliQ water
<b>Application</b>	ELISA, WB
<b>Usage recommendations</b>	ELISA- 1µg/ml, WB-1µg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

### WB image of IL 6 Receptor Monoclonal Antibody (MAB1F3)

