



Anti BSA Monoclonal Antibody

Product Details

Catalog no.	BAS-MAB02
Quantity	100µg
Clone	B11
Clonality	Monoclonal
Conjugate/ Tag	Unconjugated
Purity	Purified by protein A affinity chromatography
Form	Lyophilized/ Liquid
Buffer	Phosphate buffered saline (pH 7.4)
Immunogen	Natural Bovine serum Albumin
Applications	ELISA, WB
Usage recommendations	ELISA- 1µg/ml, WB-1µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

WB Image of Anti BSA Monoclonal Antibody

